KULLANMA TALİMATI

Optive FUSION® UD Göz Damlası

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 AbbVie Logistics B.V.,
Zuiderzeelaan 53, Zwolle,
8017 JV, Hollanda

 devices.abbvie.com

 ALLERGAN PHARMACEUTICALS
INTERNATIONAL LIMITED
Clonshaugh Business
& Technology Park,
Dublin 17, D17 E400, İrlanda

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20077513

abbvie		Date: 28 Feb 2023	Time: 16:49	Print Colors	Non-Print Colors
Component Code	20077513	Proof	01	Black	Keyline
Site	Allergan Westport	Drawing N°	0361201		
Smallest Body Text Size	9 pt	Braille code	N/A		
Site Packaging Barcode N°	N/A				
Notes	PDF proof output to A3 laser at 85%.				

INSTRUCTIONS FOR USE

Optive FUSION™ UD EYE DROP

Intended Purpose and Indications

OPTIVE FUSION™ UD Eye Drop is intended to lubricate and hydrate the surface of the eye(s) and is indicated to help provide temporary relief of dryness, irritation, burning, discomfort and/or any other symptoms caused by dry eye disease (DED) or environmental factors.

Contents of the packaging

OPTIVE FUSION™ UD Eye Drop is packaged within a plastic vial, flow wrapped (foil pouched), with an insert, and then in a card carton. OPTIVE FUSION™ UD Eye Drop solution contains sodium carboxymethylcellulose 0.5%, glycerin 1.0%, sodium hyaluronate 0.1%, levocarnitine, erythritol, sodium lactate, potassium chloride, calcium chloride dihydrate, magnesium chloride hexahydrate, purified water and hydrochloric acid or sodium hydroxide for pH adjustment. Product does not contain any preservatives.

Clinical Benefits

OPTIVE FUSION™ UD Eye Drop provides instant soothing relief of the feeling of dryness, irritation, burning, discomfort and/or any other symptoms caused by dry eye disease (DED) or environmental factors.

Performance Characteristics

OPTIVE FUSION™ UD Eye Drop can:

- Reduce eye discomfort, eye dryness and eye irritation
- Improve clinical signs of dry eye disease and tear film stability

OPTIVE FUSION™ UD Eye Drop is safe for use on the ocular surface in general and can be used repeatedly as needed, without any minimal time interval between two instillations.

OPTIVE FUSION™ UD Eye Drop is compatible with all types of lens materials (soft and rigid gas-permeable contact lenses).

OPTIVE FUSION™ UD Eye Drop is free of preservatives and well tolerated by patients with sensitive eyes and can be used by post-operative users (e.g. following LASIK surgery).

Target Patient Population and Intended Users

OPTIVE FUSION™ UD Eye Drop is an over-the-counter (OTC) ophthalmic product that is intended to be used by lay persons who

may experience dryness, irritation, burning, discomfort and/or any other dry eye symptoms and does not require specific training for application.

The device has no restrictions on the target treated population, except the one listed in the Contraindications. However, safety and effectiveness have not been demonstrated in pediatric patients. There are no specific study data on the use of the products during pregnancy and lactation in humans.

Contraindications

OPTIVE FUSION™ UD Eye Drop is contraindicated in patients with hypersensitivity to any ingredients in this product.

Warnings and Precautions

- Do not swallow solution.
- To avoid contamination or possible eye injury, do not touch vial tip on any surface, and avoid direct contact with the eye.
- Keep out of the reach of children.
- Do not use if solution changes color or becomes cloudy.
- Do not use after expiry date marked on the product.
- For single use only.
- Do not use if the unit dose vial is damaged
- Store at room temperature
- Allow 5 minutes between the administration of other ophthalmic products.
- Keep unused product in pouch.
- This product is not intended to be re-used and should be discarded immediately after initial use.
- This product is preservative free and should be used immediately after opening to avoid contamination which could result in eye infections.
- Discontinue use and consult a doctor if you experience eye pain, changes in vision, continued redness or irritation of the eye, or if the condition worsens.
- May cause transient blurring of vision which may impair the ability to drive or operate machines. The patient should wait until their vision has cleared before driving or using machinery.

Undesirable Side-Effects

OPTIVE FUSION™ UD Eye Drop clinical studies and post-marketing experience have identified the following residual risks:

- Eye irritation, Conjunctival hyperaemia, Eye redness, Foreign body sensation in eye, Erythema of eyelid, Blepharitis, Abnormal sensation in eye, Eye pain, Vision blurred, Punctate keratitis, Diffuse lamellar keratitis, Eye pruritus
- Superficial injury of the eye

Directions

- Ensure the vial is intact
- To open completely twist off the tab
- Place 1 or 2 drops in the affected eye(s) as needed
- Store unused vials in the original carton
- Store unused vials in the pouch
- Discard any remaining solution, and dispose of the container properly after use
- If used for post-operative care (e.g., following LASIK surgery), it is recommended for the patient to follow their eye doctor's instructions.
- Any unused product or waste material should be disposed of in accordance with local regulation. No special requirements for disposal.

Reporting

Report any suspected product related adverse reaction to Safety.Turkey@abbvie.com and report product related events to Turkish Medicines and Medical Devices Agency.

How supplied

OPTIVE FUSION™ UD Eye Drop is supplied in 30 pack of 0.4 mL unit dose vials, and each carton contains vials of the product, flow wrapped (foil pouched), and an insert with instructions for use.

ALLERGAN SALES, LLC,
2525 Dupont Drive,
IRVINE CA 92612,
USA

AbbVie Logistics B.V.,
Zuiderzeelaan 53, Zwolle,
8017 JV, Netherlands

devices.abbvie.com

0344

EC REP

ALLERGAN PHARMACEUTICALS
INTERNATIONAL LIMITED
Clonshaugh Business
& Technology Park,
Dublin 17, D17 E400, Ireland

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abbvie

© 2023 Allergan, Inc.
™ & ® marks owned By Allergan, Inc.

abbvie		Date: 28 Feb 2023	Time: 16:49	Print Colors	Non-Print Colors
Component Code	20077513	Proof	01	Black	Keyline
Site	Allergan Westport	Drawing N°	0361201		
Smallest Body Text Size	9 pt	Braille code	N/A		
Site Packaging Barcode N°	N/A				
Notes	PDF proof output to A3 laser at 85%.				

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AbbVie Tibbi İlaçlar San. ve Tic. Ltd. Şti.